

Finally after 50 Years...



STEREOPSIS TESTS

The FLY Stereo Acuity Test

NOW with **lea** SYMBOLS®

Finally after
+50 Years...
NEW !!!
FLY TEST



P/N 1000

The Butterfly Stereo Acuity Test

NOW with **lea** SYMBOLS®

Finally after
25 Years...
NEW !!!
BUTTERFLY TEST



P/N 1010

Tests designed to rapidly test for amblyopia & strabismus using gross thru fine stereopsis
(4800 to 20 sec of arc)

- ❖ GRADED CIRCLE TEST FROM 400 SEC **NOW DOWN TO 20 sec**
- ❖ NEW TECHNOLOGY → NO MONOCULAR CLUES
- ❖ NEW IMPROVED BOOKLETS
- ❖ ANSWER KEY ON BACK COVER
- ❖ **NOW INCLUDES** ADULT & PEDIATRIC POLARIZED VIEWERS

Random Dot Stereo Acuity Test

NOW with **lea** SYMBOLS®

Finally after
30 Years
NEW !!!
RANDOM DOT
TEST



P/N 1005

Random Dot E Depth Perception Test

Finally after
+30 Years...
NEW !!!
RANDOM DOT E
TEST



P/N 1015


Designed to rapidly test for amblyopia & strabismus in early & non-readers and non-verbal children & adults

- ❖ EXPANDED RANDOM DOT LEA TEST
(500, 250, 125, 63 sec of arc)
- ❖ GRADED CIRCLE TEST **NOW DOWN TO 12.5 sec**
- ❖ NEW TECHNOLOGY → NO MONOCULAR CLUES
- ❖ NEW IMPROVED BOOKLETS
- ❖ ANSWER KEY ON BACK COVER
- ❖ **NOW INCLUDES** ADULT & PEDIATRIC POLARIZED VIEWERS
- ❖ CHOICE OF DOT SIZE ACUITY
- ❖ NEW IMPROVED TEST FRAMES
- ❖ INCLUDES ADULT & **NEW KID-FRIENDLY PEDIATRIC VIEWERS**



PROJECTOR SLIDES

Designed for use in all manual Chart Projectors including products from:
R.H. Burton, Marco, Reichert, Topcon and Woodlyn



1077 PEDIATRIC ACUITY

- ❖ Allen Acuity 20/200 - 20/30 (6 Lines)
- ❖ Tumbling E 20/200 - 20/20 (11 Lines)
- ❖ Landolt C 20/200 - 20/20 (8 Lines)

1103 VERGENCE-ACUITY

- ❖ Tumbling E 20/400 - 20/200 (3 Lines)
- ❖ Sloan Letters 20/100 - 20/10 (14 Lines)
- ❖ Astigmatic Clock
- ❖ Vertical Lines for Vergence 20/40 - 20/10 (6 Lines)
- ❖ Sloan Letters 20/50 - 20/10 (6 Lines)
- ❖ Astigmatic Blocks

1113 ASTIGMATISM-ACUITY

- ❖ Sloan Letters 20/300 - 20/10 (14 Lines)
- ❖ Astigmatic Clock
- ❖ Sloan Numbers 20/50 - 20/10 (6 Lines)
- ❖ Tumbling E 20/200 - 20/15 (7 Lines)

1180 ADULT ACUITY

- ❖ Sloan Letters 20/400 - 20/10 (16 Lines)
- ❖ Sloan Numbers 20/70 - 20/10 (8 Lines)
- ❖ Tumbling E 20/200 - 20/15 (11 Lines)

1192 MULTI ACUITY

- ❖ Sloan Letters 20/400 - 20/15 (15 Lines)
- ❖ Tumbling E 20/200 - 20/25 (9 Lines)
- ❖ Allen Acuity 20/200 - 20/30 (6 Lines)

1193 FAMILY PRACTICE

- ❖ Tumbling E 20/400 - 20/200 (3 Lines)
- ❖ Sloan Letters 20/100 - 20/10 (11 Lines)
- ❖ Astigmatic Clock
- ❖ Sloan Letters 20/50 - 20/10 (6 Lines)
- ❖ Tumbling E 20/50 - 20/15 (5 Lines)
- ❖ Allen Acuity 20/100 - 20/40 (4 Lines)

3200 LEA SYMBOLS® & LETTERS

- ❖ LEA Symbols® 20/125 - 20/10 (19 Lines)
- ❖ (Multiple Lines 20/50 - 20/10)
- ❖ Sloan Letters 20/400 - 20/10 (20 Lines)
- ❖ (Multiple Lines 20/50 - 20/10)

3201 LEA NUMBERS® & LETTERS

- ❖ LEA Numbers® 20/125 - 20/10 (17 Lines)
- ❖ (Multiple Lines 20/50 - 20/10)
- ❖ Sloan Letters 20/400 - 20/10 (20 Lines)
- ❖ (Multiple Lines 20/50 - 20/10)

3202 LEA SYMBOLS/NUMBERS® & LETTERS

- ❖ LEA Symbols® 20/100 - 20/10 (18 Lines)
- ❖ LEA Numbers® (16 Lines)
- ❖ (Multiple Lines 20/50 - 20/10)
- ❖ Sloan Numbers 20/100 - 20/10 (16 Lines)
- ❖ (Multiple Lines 20/50 - 20/10)

3203 LEA SYMBOLS® & LETTERS

- ❖ LEA Symbols® 20/200 - 20/10 (20 Lines)
- ❖ (Multiple Lines 20/50 - 20/10)
- ❖ Sloan Letters 20/250 - 20/10 (20 Lines)
- ❖ (Multiple Lines 20/50 - 20/10)

3204 LEA SYMBOLS® & LETTERS

- ❖ LEA Symbols® 20/200 - 20/16 (11 Lines)
- ❖ Sloan Letters 20/400 - 20/16 (17 Lines)
- ❖ (Multiple Lines 20/40 - 20/20)

3205 LEA SYMBOLS/NUMBERS® & LETTERS

- ❖ LEA Symbols® 20/125 - 20/16 (10 Lines)
- ❖ LEA Numbers® 20/125 - 20/16 (10 Lines)
- ❖ Sloan Letters 20/125 - 20/16 (14 Lines)
- ❖ (Multiple Lines 20/40 - 20/20)

3206 LEA NUMBERS® & LETTERS

- ❖ LEA Numbers® 20/200 - 20/16 (11 Lines)
- ❖ Sloan Letters 20/125 - 20/16 (17 Lines)
- ❖ (Multiple Lines 20/40 - 20/20)

3207 LEA SYMBOLS® & LETTERS

- ❖ LEA Symbols® 20/200 - 20/16 (11 Lines)
- ❖ Sloan Letters 20/400 - 20/16 (17 Lines)
- ❖ (Multiple Lines 20/40 - 20/20)

1234 VECTOGRAPHIC ADULT

- ❖ Sloan Letters Monocular 20/100 - 20/15 (9 Lines)
- ❖ Astigmatic Clock
- ❖ Acuity Balance 20/30 - 20/20 (3 Lines)
- ❖ Acuity Suppression 20/40 - 20/25 (3 Lines)
- ❖ Sloan Letters Binocular 20/30 - 20/15 (4 Lines)
- ❖ Fixation Disparity (with/without Fusion Lock)
- ❖ Stereopsis 260 - 73 seconds of arc

1246 VECTOGRAPHIC PEDIATRIC

- ❖ Training Demonstration
- ❖ Tumbling E Monocular 20/100 - 20/15 (9 Lines)
- ❖ Acuity Balance 20/30 - 20/20 (3 Lines)
- ❖ Tumbling E Binocular 20/30 - 20/15 (4 Lines)
- ❖ Fixation Disparity with Fusion Lock
- ❖ Stereopsis 215 - 40 seconds of arc

