


Contact Lens Tools

GuldenOphthalmics



Contact Lens Check & Keratometer Calibrator

- ☒ Easily Mounts on Most Keratometers & Slit Lamps
- ☒ Check for Warpage or Any Changes to the Base Curve & Calibrate your Keratometer
- ☒ Eliminates a need for a Radiuscope
- ☒ Economical 15179 \$49



Keratometer Range Extender

- ☒ Extends range; flatter to 30.87D or steeper to 60.62D.
- ☒ Attaches easily for temporary or permanent use
- ☒ Eliminates, bothersome taping & alignment of lenses
- 15180 \$39



Contact Lens Polisher Kit

- ☒ Cordless
- ☒ Palm-size
- ☒ Fast motor saves time
- ☒ Includes flat and cone shaped foam pads
- 15172 \$99

225 Cadwalader Avenue
Elkins Park, PA 19027-2020
215-884-8105

Fax 215-884-0418

e-mail: info@guldenophthalmics.com

www.guldenophthalmics.com

12/27/2011