



**Gulden Ophthalmics**

225 Cadwalader Avenue

Elkins Park, PA 19027-2020, USA

Web: [www.guldenophthalmics.com](http://www.guldenophthalmics.com)

Email: [info@guldenophthalmics.com](mailto:info@guldenophthalmics.com)

Phone: +1(215) 884-8105



**Qarad B.V.**

Flight Forum 40

5657 DB Eindhoven

The Netherlands



---

**Instructions:**

1. To raise or lower, **ALWAYS** hold the base with one hand and turn the inner tube with the other hand.
2. Raise the inner tube to desired height twist clockwise to lock.
3. Twist counter clockwise to unlock and readjust.