

Gulden Ophthalmics

225 Cadwalader Avenue, Elkins Park, PA 19027
(800) 659-2250 Fax (215) 884-0418 e-mail: info@guldenophthalmics.com
www.guldenophthalmics.com
Time Saving Tools

This document is to certify the following:

Orbital Implants (Eye Spheres) : 14150, 14151, 14152, 14153, 14154, 14155, 14156,
14250, 14251, 14252, 14253, 14254, 14255, 14256

Conformers : 14164, 14165, 14166, 14167, 14168, 14169,
14264, 14265, 14266, 14267, 14268, 14269

Orbital Implants and Conformers are made of polymethylmethacrylate (PMMA) acrylic and are Magnetic Resonance Imaging (MRI) safe and compatible.

These implants are completely non-magnetic, non-electrically conductive, and non-RF reactive.

Appropriate records are available in company files.

 5/10/2016

Authorized Signature & Date

Steve McLaughlin

Q.C. Manager

